



# ACCOUNT GENERAL INFORMATION & PREFERENCES

Thank you for choosing to work with Vitality Dental Arts. We have learned that all doctors have varying degrees of preferences when it comes to how their restorations are fabricated. Learning these preferences allows us to deliver consistent and predictable results. Please take a minute to fill out the information below. Thank you!

**Doctor Name:** \_\_\_\_\_

**Address:** \_\_\_\_\_

**Phone #:** \_\_\_\_\_ **E-Mail:** \_\_\_\_\_

**How did you hear about us:**  
\_\_\_\_\_

**Office Hours:** MON \_\_\_\_\_ TUE \_\_\_\_\_ WED \_\_\_\_\_ THU \_\_\_\_\_ FRI \_\_\_\_\_

**Preferences:** These guidelines will be your defaults for every case sent to us unless otherwise indicated on the prescription form. You can change these preferences at anytime you wish by calling in or requesting a new form.

**Contacts:** (our default is broad & normal. We verify all contacts on a solid model)

- Default  Broad  Point  Normal  Tight  Light

**Occlusal Contact:** (our default is 0.5mm out of occlusion)

- Default  1mm out of occlusion  Positive contact  Other: \_\_\_\_\_

**Tight Bites:** (many doctors allow us to make minor adjustments when needed – choose your preference.)

- Ok to adjust and mark opposing if less than 0.5mm clearance  Ok to adjust and mark the prep if less than 0.5mm clearance  Make a reduction coping  Send back for a new impression  Always call

**Undercut or Unparallel Preps:** (many doctors allow us to make minor adjustments when needed – choose your preference.)

- Ok to adjust and mark prep(s) to allow the restoration to draw  Send back for a new impression  Always call

**Preferred PFM Collar Design:**

- Standard 1mm lingual  Micro 0.5mm lingual collar  No collar  Other: \_\_\_\_\_

**Preferred PFM Metal Type:**

- Non-precious  Semi- Precious  High Noble

**Occlusal Staining:**

- None  Light  Medium  Dark

**Are there any other preferences that you would like to bring to our attention:**

\_\_\_\_\_  
\_\_\_\_\_  
\_\_\_\_\_