



# IMPLANTPACK™

POWERED BY \* VITALITY DENTAL ARTS

IMPLANT 1-800-399-0705  
FIXED Rx www.implantpack.com

Dr. Name \_\_\_\_\_

Address \_\_\_\_\_

City \_\_\_\_\_ State \_\_\_\_\_ Zip \_\_\_\_\_

Phone # \_\_\_\_\_ Email \_\_\_\_\_

Patient Name \_\_\_\_\_

DATE DUE \_\_\_\_\_


BY 5:00 PM

LAB USE

## IMPLANT SYSTEM:

\_\_\_\_\_

### INDICATE IMPLANT DIAMETER BELOW:



**Rx** CASE INSTRUCTIONS  
(additional space on back)

Date \_\_\_\_\_

Signature \_\_\_\_\_ License # \_\_\_\_\_


## SELECT YOUR IMPLANTPACK™

☐ VITALITY IMPLANT PACK (FDA, 510K CLEARED) ☐ ATLANTIS IMPLANT PACK ☐ 3I ENCODE EMPOWERED IMPLANT PACK ☐ AUTHENTIC BRAND IMPLANT PACK

☐ CEMENT-RETAINED OR ☐ SCREW-RETAINED

### SELECT RESTORATION

☐ Vitality Solid Zirconia ☐ Vitality Esthetic ☐ 3I ENCODE EMPOWERED IMPLANT PACK ☐ AUTHENTIC BRAND IMPLANT PACK  
☐ Vitality Layered Zirconia ☐ IPS e.max ☐ Crown ☐ Splinted Crowns ☐ Bridge ☐ Vitality Solid Zirconia (w/ Ti-Base)  
☐ Non-Precious PFM ☐ Noble PFM ☐ High Noble PFM ☐ Vitality Esthetic (w/ Ti-Base)  
☐ Titanium ☐ Gold-Colored Titanium ☐ Prepare stock abutment ☐ Vitality Layered Zirconia  
☐ Zirconia ☐ IPS e.max (w/ Ti-Base)  
☐ Non-Precious PFM  
☐ Noble PFM  
☐ High Noble PFM



### SHADE INFORMATION


Final Shade \_\_\_\_\_

Occlusal Staining: ☐ Light\* ☐ Med ☐ Dark ☐ None

### PARALLEL ABUTMENTS

☐ No ☐ Yes (indicate which abutments will have restorations splinted together for insertion)

### ABUTMENT MARGIN DEPTH




If left blank, default values will be used

### ABUTMENT MARGIN DESIGN



### ABUTMENT EMERGENCE PROFILE



### CONTOUR AND OCCLUSION DESIGN

Embrasures: ☐ Closed\* ☐ Open  
 Occlusion: ☐ Light\* ☐ Ideal ☐ Open \_\_\_\_ mm ☐ Out  
 Contacts: ☐ Broad & Tight\* ☐ Pinpoint ☐ Light

### RESTORATION PONTIC DESIGN



### CERAMIC METAL DESIGN



\*Standard unless specified otherwise

### IF NO OCCLUSAL CLEARANCE

☐ Call doctor ☐ Spot opposing ☐ Metal occlusion  
☐ Metal island ☐ Make this a permanent note in my master file

# SENDING CASES

## TRADITIONAL IMPRESSIONS

1. Call 1-800-399-0705 or log in to your online portal at [www.vitalitydentalarts.com](http://www.vitalitydentalarts.com) to request shipping supplies, including prescriptions, boxes, and preprinted shipping labels.



2. Separate the impressions, models, and Rx into individual bags for accurate identification. You can combine multiple cases in one box. Seal the box securely with packing tape and attach the shipping label.



3. Call 1-800-399-0705 or log in to your online portal at [www.vitalitydentalarts.com](http://www.vitalitydentalarts.com) to request an in-office pickup. We will send one of our drivers, a local courier service, or UPS based on your location.



## DIGITAL IMPRESSIONS



iTero



3shape TRIOS



SIRONA



Lava 3M ESPRIMO  
Oralmax One Scanner C.B.S.



Carestream  
DENTAL



PLANMECA  
E4D TECHNOLOGIES



dental wings



CONDOR  
DENTAL

We work with ALL major impression scanners: Cadent iTero, 3Shape Trios, Sirona Cerec, 3M True Definition/Lava COS, Carestream, E4D Planmeca, Dental Wings, and Condor. We are one of the nation's top-producing labs for processing digital impressions.

No other laboratory has put as much effort and resources into understanding and preparing their workflow for the upcoming digital revolution as we have. Our dental laboratory's foundation was built upon CAD/CAM Dentistry, so why not work with the best?

In addition to accepting all digitally scanned files, we are also an authorized distributor for all 3Shape Trios Scanners. Call today to find out how to demo a unit!